

A GUIDE TO
AVIZOR


CONTACT LENS CARE SOLUTIONS
& OCULAR HEALTH PRODUCTS



SUCCESSFULLY MANUFACTURING CONTACT LENS CARE SOLUTIONS AND OCULAR HEALTH PRODUCTS SINCE 1981

Avizor is a global leader in supplying products for the cleaning and maintenance of all types of contact lenses along with an ocular health range. Launched in 1981 as a family owned business with its headquarters and production centre in Madrid, Spain, the company started with the production of private label products but soon launched its own Avizor brand.

Currently the company exports close to 90% of its production to more than 80 countries, and employs more than 200 people around the world.



United Kingdom

EUROPE: Spain, Portugal, France, Germany, Russia, Holland, Austria, Poland

LATAM: Brazil, Mexico, Argentina, Colombia, Costa Rica

ASIA: China, Indonesia, Malaysia, Philippines, Taiwan, Vietnam

MENA: Saudi Arabia, Iran, Algeria, Turkey, Morocco, Jordan

OUR PRODUCTS

PEROXIDES

ALL IN ONE

GP'S

OCULAR HEALTH DROPS

SALINE

OCULAR HEALTH OTHERS



OUR PORTFOLIO

Avizor has one of the most extensive ranges in the market of contact lens care products and eye care solutions. Avizor continues to develop new products that meet the needs of a demanding and changing market.



+80
COUNTRIES



200
EMPLOYEES
around the world



+11000m²
50M bottles
GMP licence
Versatile logistics



Brand
Avizor brand
Private labels
Premium Products



Own
Laboratory
QA
Regulatory Affairs
R+D



Original
Own moulds/
own lens cases
Flexibility in formulas
and sizes



Complete Range
Contact lens Care
Eye Care

PEROXIDE



THE ULTIMATE INNOVATION IN CONTACT LENS CARE

PRESERVATIVE
FREE



- One step peroxide system.
- Ensures lipid cleaning thanks to a new biotechnological agent.
- Deep removal of proteins thanks to a proteolytic enzyme.
- Guaranteed disinfection against bacteria, fungi.
- Disinfecting efficacy against Acanthameoba castellanni and polyphaga.
- With Povidone, a lubricant that improves comfort to the user.
- Safe colour indicators show the tablet has been added.
- For all contact lens materials; Soft (silicone hydrogel) and RGP (Ortho-k).



ACTION
IN 2 HOURS



COMPLETE
REMOVAL
OF PROTEINS



NO NEED
TO RUB

Peroxide (%): 3% • **Tablet:** Biotechnological agents (lipid and protein cleaning), Povidone (lubricating agent), catalase (neutralising) and colour indicators.

Formats: 2 x 350 ml (90 days) / 350 ml (45 days) / 225 ml (30 days) / 60 ml (8 days)



PEROXIDE



THE MOST INNOVATIVE CARE SYSTEM FOR SCLERAL LENSES

PRESERVATIVE
FREE



- One step peroxide system with a new specific lenscase lens case for scleral lenses with large diameter 23 mm and sagittal height 9 mm.
- Guaranteed disinfection against bacteria, fungi.
- **Disinfecting efficacy against Acanthamoeba castellanii and polyphaga.**
- Deep removal of proteins thanks to a proteolytic enzyme.
- With Povidone, a lubricant that improves comfort to the user.
- Safe colour indicators show the tablet has been added.
- For all contact lens materials; Soft (silicone hydrogel) and RGP (Ortho-k).



ACTION
IN 2 HOURS



COMPLETE
REMOVAL
OF PROTEINS



NO NEED
TO RUB

Peroxide (%): 3% • **Tablet:** Proteolytic enzyme, Povidone (lubricating agent), catalase (neutralising) & colour indicators.

Formats: 300 ml (30 days) / 80 ml (6 days)

ALL
IN ONE



alvera®

A UNIQUE CONTACT LENS SOLUTION CONTAINING ALOE VERA



- Specially designed for Silicone Hydrogel lenses
- Enhanced lipid removal action
- Improves cell viability
- Prevents infection and enhances UV protection.



PROTECTION



SILICONE
HYDROGEL



EFFICIENT
AGAINST
LIPIDS

Preservative (%): Polyhexamethylene biguanide 0.0002 • **Active Agent:** Poloxamer / Aloe Vera • **Lubricating agent:** Polyvinylpyrrolidone

Formats: 3 x 240 ml / 350 ml / 100 ml / 60 ml • Available for delivery direct to patient via Adaro.

ALL
IN ONE



BENEFITS OF ALVERA

Multipurpose solution for all types of soft contact lenses, and specially developed for use with silicone-hydrogel lenses.

FUNCTION	COMPONENTS	USER BENEFITS
Patented Cellular Epithelium protection.	Aloe Vera	Aloe Vera prevents bacterial infections, and protects against UV rays and eye damage.
Reduces lipid affinity that adhere to the lens surface.	Aloe Vera	Increased comfort due to less deposits.
Lubricates contact lenses and the eye.	PVP, Aloe Vera	Potent lubricant that provides all-day comfort in a natural way.
Removes deposits and debris that adhere to the lens (including lipids).	Cleaner and surfactant, Poloxamer	Increases wearing time and visual acuity while keeping your lenses cleaner.*
Efficient Disinfection and Cleaning. Conserves lenses when not in use.	Polyhexanide	One of the least toxic products without losing its efficiency during disinfection. End user safety, less risk of development of bacterial and ocular infections.

CLINICAL PAPERS

* Evaluation of Safety and efficacy of a new multipurpose disinfecting solution on SiHy CL. <http://www.journalofoptometry.org/>
Cytotoxic Effect on Corneal Surface. <https://www.omicsonline.org/>

ALL
IN ONE



Unica[®] sensitive



PRESERVATIVES



WETTING



Avizor Unica Sensitive has been especially formulated to bring extra comfort. The combination of an extraordinary lubricating agent (high molecular weight sodium hyaluronate) and the minimum concentration of preservatives, make Unica Sensitive the best option for users of soft contact lenses with sensitive eyes.

Suitable for all contact lens modalities: can be used for breaks throughout the day in which the user removes their lenses whilst: resting, in the pool, playing sports.

Some of the advantages of sodium hyaluronate are the following:

1. Vision improvement: reduces the symptoms of dry eye.
2. Increases tear stability as it reduces the evaporation of the tear film.
3. Increases hours of satisfactory use thanks to its viscoelastic properties.

Preservative (%): Polyhexamethylene biguanide 0.0001 • **Active Agent:** Poloxamer • **Lubricating agent:** Polyvinyl-pyrrolidone / Sodium hyaluronate

Formats: 3 x 240 ml / 350 ml / 60 ml / 10 x 10 ml • Available for delivery direct to patient via Adaro.

ALL clean[®] soft

ALL
IN ONE



DISINFECTION



CLEANING



LUBRICATION

THE MOST EFFECTIVE ALL IN ONE SOLUTION

Studies carried out show that All Clean Soft is the most effective all in one solution for all users of soft contact lenses including silicone hydrogel.

1. Efficient in disinfection.
2. Efficient against Acanthamoeba: in only 4 hours it eliminates cysts and trophocytes of Acanthamoeba castellanii and polyphaga.
3. Efficient in cleaning: the Poloxamer included in the formula catches the deposits facilitating the cleaning on the surface of the lens.

Preservative (%): Polyhexamethylene biguanide 0.0002 • **Active Agent:** Poloxamer • **Lubricating agent:** Polyvinylpyrrolidone

Formats: 3 x 240 ml / 350 ml / 100 ml / 60 ml • Available for delivery direct to patient via Adaro.

GP'S



GP SYSTEM: GP CLEANER + GP CONDITIONER



cleaner



conditioner

Cleaner specifically designed for rigid and gas-permeable contact lenses.

It is a very effective anionic surfactant. Contains particles that eliminate lipid deposits on rigid and gas-permeable lenses.

Wetting, storing, disinfecting and rinsing solution for rigid and gas permeable contact lenses.

Once the lenses have been cleaned, Avizor GP Conditioner covers them with a lubricant film insuring maximum comfort.

Active Agent: Sodium Lauryl sulphate / Propylene glycol

Formats: 30 ml

Lubricating Agents: Polyvinylpyrrolidone / Hypromellose

Formats: 120 ml



Multi-purpose solution to clean and disinfect rigid and gas permeable lenses.

Multi-purpose solution for the most demanding users and those who need a higher level of cleanliness.

TECHNICAL DATA

- **Sodium chloride:** which offers the necessary osmolarity.
- **Borax/Boric:** which balances the solution's pH so that it adapts to that of the eye.
- **Poloxamer:** is the surfactant agent that cleans the lens.
- **Polyhexanide:** antimicrobial disinfectant.
- **PVP & Hypromellose:** lubricating agents.
- **EDTA:** chelator, which aids the cleaning process.

Preservative (%): Polyhexamethylene biguanide 0.0002 • **Active Agent:** Poloxamer • **Lubricating agent:** Polyvinylpyrrolidone / Hypromellose

Formats: 3 x 240 ml / 240 ml / 60 ml



lacrifresh®

RANGE OF EYE DROPS



cleaning



PRESERVATIVE
FREE



aloe



PRESERVATIVE
FREE



moisture
0,10%



CONTACT LENS CLEANER

EYE IRRITATION

DRY EYE

ENVIRONMENTAL DRY EYE

LUBRICATING / WETTING AGENT

Povidone

Aloe Vera
& Hypromellose

Sodium Hyaluronate
0.10%

OSMOLARITY (mOsm/L)

300

300

300

VISCOSITY (cP)

3-5

3-5

3-5

PRESERVATIVE

Yes

None

None

BENEFITS

- Blinking action cleans the lens on the eye (no need to remove contact lenses).
- Removes the deposits that accumulate on the lens throughout the day.
- Cleans lipids.

- Relieves the symptoms of eye irritations.
- Provides a high level of lubrication and wetting action.
- Phosphate free.
- Easy application.
- Practical format.
- OSD System.

- Relieves occasional ocular dryness due to environmental situations.
- Moisturises and lubricates the eye.
- Isotonic.
- Unidose vials and OSD System bottle options.

PACK SIZES

15 ml

10 ml

10 ml / 20 x 0,4 ml

CL FRIENDLY



lacrifresh®

RANGE OF EYE DROPS



DRY EYE

MILD TO MODERATE DRY EYE¹

MODERATE TO SEVERE DRY EYE^{2,3}

LUBRICATING / WETTING AGENT	Sodium Hyaluronate 0.20% & Glycerine & Povidone	Sodium Hyaluronate 0.30% & Glycerine & Povidone
OSMOLARITY (mOsm/L)	220	160-190
VISCOSITY (cP)	16-19	48-52
PRESERVATIVE	None	None
BENEFITS	<ul style="list-style-type: none"> Relieves the symptoms of mild to moderate ocular dryness. Hypotonic. Unidose vials and OSD System bottle options. 	<ul style="list-style-type: none"> Relieves the symptoms of moderate to severe ocular dryness. Hypotonic (more than 0,20%). Unidose vials and OSD System bottle options. Overnight use (10 ml).
PACK SIZES	10 ml / 20 x 0,4 ml	10 ml / 20 x 0,4 ml
CL FRIENDLY	✓	✓

DROPS



CLINICAL PAPERS

¹ References

OCU-DRY 0,2%:
Efficacy and Safety of 0.2% Hyaluronic Acid in the Management of Dry Eye.
<http://journals.lww.com>.
Search: 10.1097/ICL.0000000000000236.

² References

OCU-DRY 0,3%:
Safety and Efficacy of an artificial Tear Containing 0.3% Hyaluronic Acid in the Management of Moderate-to-Severe Dry Eye Disease" Eye Contact Lens. 2017 Nov;43(6): 383-388. doi: 10.1097/ICL.0000000000000284

³ References

OCU-DRY 0,3%:
Carracedo G., et al., Comparison Between Viscous Teardrops and Saline Solution to Fill Orthokeratology Contact Lenses Before Overnight Wear. Eye & Contact Lens, 2017.

SALINE



saline

PHOSPHATE FREE

Rinses, hydrates and preserves all contact lens materials (including dailies & scleral lenses).

It can also be used for enzymatic cleaning.



PRESERVATIVE
FREE



Preservative (%): Polyhexanide 0.0001

Formats: 350 ml / 30 x 5 ml

enzyme



- Effervescent enzymatic tablets.
- Eliminate all protein deposits in all types of contact lenses.
- With subtilisine A.



Active Agent: Subtilisine A (2 mg / tablet)

Formats: 10 Tablets

lipid clean



Lipid Clean SiHy, is designed to be the ideal complement for contact lens users who have problems with the accumulation of lipids and other deposits on their lenses.



Preservative (%): Polyhexanide

Formats: 60 ml

OTHERS



AVIZOR

**FOR MORE INFORMATION CONTACT
YOUR AVIZOR REGIONAL SALES MANAGER
OR YOUR USUAL WHOLESALE SUPPLIER**

Avizor UK Ltd - 2 Chester Row- London SW1W 9JH
AVIZOR UK: Freephone number: 08 000 856 559
EMAIL: uk.orders@avizor.com
