

# OPTIMEL<sup>®</sup>

## Eye Care

### Utilising Manuka Honey for Chronic Dry Eyes

Developed using proprietary pharmaceutical grade honey from *Leptospermum* spp, Optimel™ eye care treatments are effective as adjunctive therapies for MGD, reducing conjunctival and lid margin redness and improving Meibomian Gland secretion quality\*.



Moderate to severe dry eye  
(Preservative free)



Mild to moderate  
dry eye



Dermatological  
conditions

Dr Keyur Patel,  
Optometrist & Clinical Director at TKS

Why do you recommend Optimel  
to your patients?

*"Optimel helps reduce the amount of inflammatory markers on the ocular surface, so the eye is less irritated, given our patients more comfortable eyes and improved tolerance to their contact lenses."*

	Optimel Drops	Optimel Gel	Optimel Cream
 <b>Concentration</b>	16.5% Manuka Honey	98% Manuka Honey	16.5% Manuka Honey
 <b>Bottle Size</b>	10ml bottle	10g tube	15g tube
 <b>Discomfort level</b>	Mild to moderate Dry Eye	Moderate to severe Dry Eye	Eczema, Dermatitis & Dry Skin
 <b>Shelf life</b>	Min. 2 years. Use within 2 months of opening	Min. 2 years. Use within 2 months of opening	Min. 2 years. Use within 3 months of opening

Best dispensed to your persistent anterior Blepharitis patients or end stage MGD patients. A great weapon for severe Dry Eye patients who are still symptomatic despite being on a high-grade PF Hyaluronic drop.